



Victorian Inclusion Agency Specialist Equipment Library Catalogue June 2018


Inclusion Support Programme (ISP) is funded by the Australian Government Department of Education and Training. As part of the ISP, the Victorian Inclusion Agency is led by Community Child Care, and is delivered in partnership with Yooralla and KU Children's Services.

Contents page



1.	Chairs and seating	Page 3
2.	Standing frames	Page 16
3.	Hoists	Page 24
4.	Toileting	Page 25
5.	Tables	Page 30
6.	Bicycles	Page 33
7.	Bead and mirrored stands	Page 34
8.	Switch and other toys	Page 38
9.	Grab bars and suction rails	Page 50
10.	Positioning systems	Page 51
11.	Ramps	Page 53
12.	Accessories	Page 55



**ALL EQUIPMENT SHOWN IN THE CATALOGUE IS A GUIDE ONLY AND IS DEPENDENT ON AVAILABILITY AT THE TIME
THE EQUIPMENT REQUEST IS MADE BY A SERVICE.**



1. Chairs and seating

Chairs and seating		
Bumbo Chair	<p>Bumbo chairs can be used on any flat ground level surface. They are a lightweight, portable plastic seat suitable for babies four months and older.</p> <p>It is shaped to support the baby when sitting. Available in only one size.</p> <p>The baby must have independent head control before using this product.</p> <p>Equipment measurements: - H 20CM W 13CM L 13CM</p>	



DES Kelly Chair	<p>The Kelly chair is suitable for a child who needs a kinder chair with extra support. It has adjustable armrests, footplate, and pommel and tray.</p> <p>Available in pre-school and primary sizes, each of which can be adjusted for growth.</p> <p>It is height-adjustable to enable a child to sit at a kinder table and participate in group activities.</p> <p>Measurements of chair: - W 33CM H 55CM L 45CM</p>	
Fixed Wooden Stool	<p>Equipment measurement: - Fixed height L 30CM W 37CM H 30CM</p>	



<p>Full Support Swings</p>	<p>The Full Support Swing seat is suitable for a child who needs a swing seat with extra support.</p> <p>It is available in three sizes (the range accommodates children from approximately 2 years to adult) and provides trunk and head support.</p> <p>The full support swing seat will fit onto standard “pigtail” swing hooks and can be adjusted in height and angle of recline.</p>	
<p>Gravity Chair</p>	<p>A foam chair which utilises gravity to hold an individual with low postural control in a semi upright position.</p> <p>The chair features raised armrests to support the user laterally. The seat cover which is incontinent resistant can be removed for cleaning.</p> <p>Standard and mobile tilting base options are available.</p> <p>Equipment measurement: - SMALL CHAIR H 55CM W 50CM L 30CM TILT WEDGE H 12CM W 51CM L 52CM WHEEL BASE H 22CM W 50CM L 52CM</p>	



<p>Jenx Corner Seat</p>	<p>The corner seat is suitable for a child who needs support to sit straight-legged on the floor. It is available in three sizes, each of which can be adjusted for growth. The corner seat can be used with the nursery table to enable access to toys. Nursery table is also available.</p> <p>Size 1: - W 50CM x H 35CM x L 55CM</p> <p>Size 2: - W 58CM x H 57CM x L 57CM</p> <p>Size 3: - W 67CM x H 65CM x L 77CM</p>	
<p>Jenx Giraffe Chair</p>	<p>The giraffe chair is suitable for a child who needs extra seating support.</p> <p>It can be adjusted to accommodate a child from approximately 12 months to 5 years, and includes a range of adjustable positioning supports.</p> <p>The giraffe chair can be used with a tray, or pushed up to a table without a tray (depending on table height) to enable participation in group activities.</p> <p>Equipment measurements: - W 30CM H 53CM L 46CM</p>	


<p>Jenx Junior Chair</p>	<p>The Jenx Junior Chair is a fully adjustable contoured modular seating system.</p> <p>It has tilt in space with a hydraulic or powered base.</p> <p>The system provides excellent postural support for a wide range of users.</p>	
<p>Juni Chair</p>	<p>A wooden framed chair with a padded contoured vinyl seat and backrest.</p> <p>The backrest features adjustable lumbar support and recline.</p> <p>The seat height and depth are adjustable.</p> <p>The armrests are made of wood and are height adjustable.</p> <p>An optional footrest is available as shown.</p> <p>Equipment measurements: - H 47CM W 58CM L 55CM (CHAIR)</p>	



<p>Leckey Advance Chair</p>	<p>The advance chair is suitable for a child who needs maximum seating support.</p> <p>It includes a range of adjustable positioning supports.</p> <p>The advance chair is height-adjustable and can be used with a tray, or pushed up to a table without a tray to enable participation in group activities.</p> <p>Equipment measurements: - H 145CM W 60CM L 75CM</p>	
<p>Leckey Easy Seat</p>	<p>The Easy Seat is suitable for a child who needs a kinder chair with extra support.</p> <p>It has adjustable armrests, footplate, and pommel, and additional side supports are available.</p> <p>Available in size 1, 2, and 3, each of which can be adjusted for growth.</p> <p>It is height-adjustable to enable a child to sit at a kinder table and participate in group activities.</p>	



<p>Leckey Everyday Seat 1</p>	<p>Designed specifically for children or young adults with moderate needs who still require pelvic stability but do not need the same advanced level of support as those with complex needs.</p>			
<p>Leckey PAL Childrens Seat</p>	<p>An adjustable chair suitable for children with mild to moderate postural support needs. Features a wooden frame, angle adjustable backrest, adjustable seat depth, height adjustable legs, footplate, armrests, and bright upholstery.</p> <p>Available in size small and medium. A range of accessories are available including lap belt, side pads, ramped cushion, and pommel.</p> <table><tr><td><p>Size 1: Green</p><p>Seat Depth 180-240mm</p><p>Seat Width 150-250mm</p><p>Backrest Height 265mm</p><p>Seat to Floor 160-240mm</p></td><td><p>Size 2: Orange</p><p>Seat Depth 220-290mm</p><p>Seat Width 170-270mm</p><p>Backrest Height 315mm</p><p>Seat to Floor 185-285mm</p></td></tr></table>	<p>Size 1: Green</p> <p>Seat Depth 180-240mm</p> <p>Seat Width 150-250mm</p> <p>Backrest Height 265mm</p> <p>Seat to Floor 160-240mm</p>	<p>Size 2: Orange</p> <p>Seat Depth 220-290mm</p> <p>Seat Width 170-270mm</p> <p>Backrest Height 315mm</p> <p>Seat to Floor 185-285mm</p>	
<p>Size 1: Green</p> <p>Seat Depth 180-240mm</p> <p>Seat Width 150-250mm</p> <p>Backrest Height 265mm</p> <p>Seat to Floor 160-240mm</p>	<p>Size 2: Orange</p> <p>Seat Depth 220-290mm</p> <p>Seat Width 170-270mm</p> <p>Backrest Height 315mm</p> <p>Seat to Floor 185-285mm</p>			


<p>Movin`Sit Inflatable Cushion</p>	<p>An inflatable wedge which serves as a seat wedge and encourages an upright posture.</p> <p>The upper surface is bumpy and provides tactile feedback while undertaking other activities.</p> <p>Available in a standard size and a junior size as shown.</p> <p>Equipment measurements: -</p> <p>SMALL: - H 7CM W 25CM L 26CM</p>	
<p>Wooden Stool</p>	<p>Height adjustable with screws with a vinyl padded seat.</p> <p>Equipment measurements: - H 36CM W 40CM L 14CM</p>	

<p>The Star X Chair</p>	<p>The Star X Chair is suitable for a child who needs extra seating support.</p> <p>It will accommodate a child from approximately 12 months to 4 years.</p> <p>It includes a range of adjustable positioning supports including thoracic, hip, and foot supports.</p> <p>The Star X Chair can be moved up to a table (depending on table height) to enable participation in group activities.</p> <p>Equipment measurement: - H 70CM W 43CM L 60CM</p>	
<p>TREC Kinder Chair</p>	<p>The TREC chair is suitable for a child who needs a kinder chair with extra support.</p> <p>It has adjustable armrests, footplate, and pommel and is available in various sizes.</p> <p>It is height-adjustable to enable a child to sit at a kinder table and participate in group activities.</p>	

<p>Tumble Forms Floor Sitter</p>	<p>The Tumble Forms Floor Sitter is suitable for a child who needs extra support for floor sitting and cannot tolerate a fully upright position.</p> <p>It is available in three sizes (the range accommodates children from approximately 6 months to primary school age), and the seat can be positioned on the base in varying degrees of recline.</p> <p>The Tumble Forms Floor Sitter can be used with an easel table, as shown to allow access to toys and activities.</p> <p>Equipment measurements: -</p> <p>Size 1 H 43CM W 43CM L 60CM including base Size 2 H 47CM W 43CM L 65CM including base Size 3 H 60CM W 43CM L 70CM including base</p>	
--	---	--

<p>Wheelie Stool</p>	<p>The wheelie stool is a mobile, height-adjustable stool designed to enable staff to work with children without bending or twisting.</p> <p>It has a wheeled base, an upholstered seat which can swivel a full 360 degrees and gas-assisted height adjustment.</p> <p>It enables staff to work at the level of the child for feeding, seated activities, and supported standing and walking.</p>	
<p>Jenx Bee Seating Unit</p>	<p>A fully adjustable, modular seating system specifically designed for young children from 4 months through to 5 years.</p> <p>The Bee seating system is an ideal introduction to supported functional seating.</p> <p>It provides an outstanding level of versatility via fixed and height adjustable base options, along with a wide variety of head, trunk, leg and foot support options.</p> <p>Size: Age range: 4 months to 5 years Min seat / max seat depth: 110-350mm Max user weight: 25kg Equipment measurements: - H 100CM W 58CM L 87CM</p>	

<p>Jenx Ladybird Floor Sitter</p>	<p>The ladybird is a prone (forward-tilted) floor sitter suitable for a child who needs support to sit straight-legged on the floor and requires a prone angle to assume correct posture.</p> <p>It can be adjusted to accommodate a child from approximately 6 months to 3 years and includes a range of adjustable positioning supports.</p> <p>The ladybird has a tray to enable access to toys.</p> <p>Equipment measurements: - H 50CM W 60CM L 90CM</p>	
<p>Scallop Seat</p>	<p>The Scallop is suitable for children who require a little extra help while floor sitting or sitting on a chair.</p> <p>This product is also great for children with learning or behavioural disorders as it provides tactile stimulation which helps to reassure the child and enable them to relax and concentrate while using the product.</p>	

<p>Bilibo Seat</p>	<p>Deceptively simple and magically attractive, Bilibo's iconic shell shape and bright colours immediately appeal to children.</p> <p>Bilibo arouses curiosity, engages the imagination and playfully trains basic motor skills and balance.</p> <p>Both boys and girls find unlimited ways to play with Bilibo – using it to rock and spin in, hide under, sit on, tote with and peek through.</p> <p>Bilibo provides hours of fun, indoors and outdoors – all year round</p>	
--------------------	--	---

2. Standing frames

Standing frames

DES Standing Frame

This prone standing frame features multiple adjustments on all components to grow with the child.



It consists of a wooden frame, padded vinyl supports, and quick release straps with buckles.


Frame incorporates a manually operated tilt adjustment with a positive locking mechanism and a solid back support to hold child in position.


Mobile frame folds flat for storage.


Equipment measurements: - H 96CM W 55CM L 73CM





<p>DES Supine Standing Frame</p>	<p>The DES Supine Standing Frame is a supine (reclined)/upright standing frame suitable for a child who needs support to stand and may not be able to tolerate a fully upright position and includes a range of adjustable positioning supports.</p> <p>It can be adjusted to enable a standing position from horizontal to fully upright.</p> <p>The DES supine standing frame has an activity tray to allow access to toys.</p> <p>Equipment measurements: - H 106CM W 40CM L 70CM</p>	
<p>Jenx Supine Stander</p>	<p>The Jenx Supine Stander is a robust stander with comfortable, adjustable accessories to provide safe standing to suit children with a wide variety of needs. While supported and upright, children can participate in a range of activities.</p> <p>Available in two sizes to suit children from 9 months onwards.</p> <p>A robust stander with a range of accessories including trays to suit a variety of needs.</p> <p>Equipment measurements: - Size 1 H 96CM W 70CM L 106CM</p>	



<p>Dynamic Standing Frame</p>	<p>The Dynamic Standing Frame is suitable for a child who needs support to stand, and is designed to allow the child some movement.</p> <p>It is available in the small and medium size and can be adjusted in height.</p> <p>The Dynamic Standing Frame can be used with a Taddy height adjustable cut out table, or it can be pushed up to an easel, water trough or high table to allow access to activities with other children.</p> <p>Equipment measurements: - H 60CM W 43CM L 60CM</p>	
<p>EasyStand Bantam Chair</p>	<p>The Bantam line of standers are the only sit to stand standers that can combine the option of supine positioning to create infinite positioning possibilities from 90-90 sitting, fully vertical standing, full supine lying and anywhere in between.</p> <p>No other standers can provide the amount of positioning that is possible with the EasyStand Bantam.</p>	

<p>Jenx Humphrey Standing Frame</p>	<p>The Humphrey is a prone (forward-tilted) upright standing frame suitable for a child who needs support to stand and may not be able to tolerate a fully upright position.</p> <p>It can be adjusted to accommodate a child from approximately 3 to 7 years, and includes a range of adjustable positioning supports. It can be adjusted to enable a standing position from horizontal to fully upright.</p> <p>The Humphrey standing frame has an activity tray that can be used either with an insert to create a flat play surface, or with a bowl for sensory play.</p>	
-------------------------------------	---	---


<p>Jenx Monkey Standing Frame</p>	<p>The Monkey is a prone (forward-tilted) upright standing frame suitable for a child who needs support to stand and may not be able to tolerate a fully upright position.</p> <p>It can be adjusted to accommodate a child from approximately 12 months to 5 years, and includes a range of adjustable positioning supports.</p> <p>It can be adjusted to enable a standing position from almost horizontal to fully upright.</p> <p>The monkey standing frame has an activity tray that can be used either with an insert to create a flat play surface, or with a bowl to allow sensory play.</p> <p>Equipment measurements: - H 72CM W 58CM L 87CM</p>	
-----------------------------------	--	---

<p>Jenx Penguin Stander</p>	<p>The Jenx Penguin Stander has been designed to provide excellent support for children in upright standing, and the range and support of adjustment options available creates a product that is suitable for children aged 18 months to 4 years</p>	
<p>Jenx Prone Standing Frame</p>	<p>The Jenx prone standing frame is a prone (forward-tilted) upright standing frame suitable for a child who needs support to stand and may not be able to tolerate a fully upright position.</p> <p>It can be adjusted to enable a standing position from upright to 40 degrees.</p> <p>The Jenx prone standing frame has an activity tray and bowl to allow access to toys.</p> <p>Equipment measurements: - H 95CM W 65CM L 95CM</p>	



<p>Leckey Paediatric Squiggles standing frame</p>	<p>The Squiggles Stander is a bright colourful three-in-one stander for younger children.</p> <p>This versatile, adjustable stander offers prone, upright, and supine standing in one product, allowing the therapist to offer standing therapy to address a child's needs and support inclusion.</p> <p>Equipment measurements: - H 110CM W 60CM L 130CM</p>	
---	---	---


<p>Jenx Multistander – Supine / prone / upright standing system</p>	<p>This standing frame is suitable for children between the ages of nine months and six years old.</p> <p>It can be configured to support supine, prone, and upright standing with two types of headrests.</p> <p>The Multistander also allows for a tray to be attached in the supine, prone and the upright position.</p> <p>Equipment measurements when closed down: - H 80CM W 58CM L 76CM</p>	
<p>Jenx Rabbit Standers</p>	<p>The Jenx Rabbit Paediatric Mobile Stander is a secure, reliable, and easy-to-adjust standing aid that supports special needs children in the vertical and prone positions.</p> <p>The tray comes with a bowl insert.</p> <p>Equipment Measurements: - H 60CM W 38CM L 55CM</p>	



3. Hoists


Hoists		
Hoist	<p>The mobile hoist is suitable for situations where manual lifting of a child is not viable, usually for a larger child who is not walking and is unable to do a standing transfer.</p> <p>It is useful for lifts across short distances, usually to and from specific positions, for example from chair to change table or chair to floor.</p> <p>Slings are available in a variety of styles depending on the amount of support required.</p> <p>The mobile hoist is powered by a rechargeable battery and is operated with a hand-held control.</p> <p>Equipment measurements: - H 50CM W 65CM L 105CM</p>	


4. Toileting

Toileting		
Electric Change Table	<p>The electric change table is suitable for a larger child who is able to do a standing transfer to and from the change table, reducing the need for lifting.</p> <p>It has an electric height-adjustment mechanism to enable it to be lowered for transfers and raised to a suitable working height for nappy changes. Cot sides or rails provide extra security.</p> <p>The electric change table is fitted with castors to enable ease of movement for storage. Various sizes stocked in Library.</p>	
Flamingo Toilet Chair	<p>The R82 Flamingo Toilet chair provides extra toileting support.</p> <p>It is available in several sizes to accommodate children from approximately 3 years to adolescent. It has a range of adjustable positioning supports, and can tilt in space.</p> <p>It can be used as a free-standing commode or positioned over the toilet.</p> <p>Equipment measurements: - H 125CM W 48CM L 76CM (stocked in Library)</p>	



<p>High backed toilet support</p>	<p>Supplied with a harness for full trunk support.</p> <p>Smaller model is intended for children up to 10 years.</p> <p>Reduced seat depth (RSD): sits 4" forward to lessen seat depth when attached to U-Shape Mounting Base.</p> <p>Dimensions of supports: height adjustable 19-24"; width at top 8" and width at bottom 12".</p>	 <p>The image shows a white, high-backed toilet support device. It features a blue harness with two vertical straps and two horizontal straps, secured with buckles. The device is mounted on a white toilet seat. The backrest is adjustable, and the device is designed to provide full trunk support for users.</p>
-----------------------------------	--	---



<p>Hoisting slings</p>	<p>Manufactured under the highest quality control standards.</p> <p>The extremely durable slings come in a range of fabrics and sizes to suit the individual child when used in conjunction with an electric powered hoist to assist educators and staff with toileting needs. Various sizes stocked in Library.</p>	
<p>Low back toilet support</p>	<p>Designed to fit any school toilet, the toilet supports strong and secure frame allows it to be easily adjusted for maximum user control and comfort, and also allows the low back toilet support to be quickly removed when not in use for easy storage.</p>	



<p>Platform step and rails</p>	<p>This multi-purpose step up platform has support handrails on both sides to give 240mm extra height.</p> <p>It is designed to step up to toilets, washbasins, tables, and other areas where a sturdy, stable step is required. It has a slip resistant surface.</p> <p>Equipment measurements: - H 55CM W 51CM L 75CM</p>	
<p>Infant Potty Chair - Adjustable</p>	<p>Contoured fibreglass mounted seat, mounted on zinc treated, epoxy coated steel frame. Three sizes available with optional chest strap available.</p> <p>Fits most toilets when fitted with adjustable legs. Adjustable legs measurements are 60 – 240 mm height to potty chair.</p>	

<p>Rifton Blue Wave Toileting System</p> <p>(Small and large chair systems)</p>	<p>The Blue Wave Toileting System is suitable for a child who needs extra support for toileting.</p> <p>It is available in two sizes to accommodate children from approximately 3 years to adolescent, and has a range of adjustable positioning supports. It can be used with or without a mobile frame.</p> <p>With the mobile frame it can be used as a free-standing commode with a potty attachment, or as an over-toilet commode without the potty attachment.</p> <p>Without the mobile frame the seat can be mounted on the toilet.</p> <p>Equipment measurements: - H 105CM W 57CM L 70CM</p>	
---	--	---


5. Tables

Tables		
DES Height Adjustable Easel Table	<p>The height adjustable easel table provides a slanted surface for activities such as drawing and reading.</p> <p>It is angle-adjustable, bringing activities into the child's line of vision, improving posture and encouraging fine motor development.</p> <p>Equipment measurements when closed down: -</p> <p>H 26CM W 62CM L 42CM</p>	
Height adjustable table	<p>Height adjustable with fold down legs to provide flexibility with ensuring the table is at the correct height for the child.</p> <p>Equipment measurements: - H Adjustable W 115CM L 70CM</p>	


<p>Height adjustable table (with wheelchair cut out)</p>	<p>Height adjustable table with left or right hand cut out for power wheelchair control.</p> <p>The Library stocks a variety of table sizes.</p>	
<p>Jenx Nursery Tables</p>	<p>The Jenx Nursery Table is available in two sizes to work with the Jenx corner chairs.</p> <p>It's a sturdy table with a large, clear surface area that is big enough for a child to play or work at most activities with ease. The table (and its surface) are ideal for drawing and writing activities.</p> <p>Equipment measurements: - H 45CM W 70CM L 50CM</p>	



<p>Four person non-adjustable table</p>	<p>Pink Non-adjustable table.</p> <p>Equipment measurements: - H 51CM W 50CM L 130CM</p>	
<p>Taddy Table</p>	<p>Height adjustable legs with table cut out provides flexibility when used in conjunction with standing frames or sitting systems.</p> <p>Equipment measurements: - H 58CM W 50CM L 50CM</p>	

6. Bicycles

Bicycles		
Supportive Tricycle	<p>The Supportive Tricycle is suitable for a child who is unable to use a regular tricycle.</p> <p>It can be adjusted to accommodate a child from approximately 3 to 6 years, and includes a range of adjustable positioning supports.</p> <p>The Library stocks a variety of sizes.</p>	

7. Bead and mirrored stands


Bead and mirrored stands		
Somatosensory Bead Chain stand with lights and music	<p>The Somatosensory Bead Chain is designed to provide hours of entertainment and tons of smiles.</p> <p>Paediatric tactile stimulation is important for a number of therapy modalities offer necessary tactile input, increasing tactile awareness as well as improving the sense of touch.</p>	


<p>Mirror with beads - small</p>	<p>Comes with a peek-a-boo mirror! Push aside the bead chain, and you'll see yourself in the unbreakable plastic mirror.</p> <p>Great for providing visual stimulation and body identification.</p>	
<p>Mirror with beads - medium</p>	<p>Comes with a peek-a-boo mirror! Push aside the bead chain, and you'll see yourself in the unbreakable plastic mirror.</p> <p>Great for providing visual stimulation and body identification.</p>	

Mirror with
beads - **large**



Comes with a peek-a-boo mirror! Push aside the bead chain, and you'll see yourself in the unbreakable plastic mirror.



Great for providing visual stimulation and body identification.







<p>Multi bead chain lights and music</p>	<p>Visual and sensory stimulation toy.</p> <p>Movement of the beads will turn on a built-in music box.</p> <p>As long as the beads are making contact with the metal bar, the music stays on.</p> <p>Size: H 32CM W 37CM L 27CM</p> <p>Weight: 1.02kg</p> <p>Requires 2 AA Batteries</p>	
--	--	---



8. Switch and other toys


Switch and other toys		
Big track ball	<p>The large ball requires less fine motor control than a standard trackball and it is ruggedly built.</p> <p>It has a left and a right mouse click button located behind the trackball to avoid unwanted mouse clicks.</p>	
5 Function Activity Centre	<p>This device increases visual attention, tactile awareness, colour recognition, and auditory development.</p>	



Casio keyboard	Musical keyboard to encourage fine motor skills and sensory interaction	
Music blocks	<p>Music Blocks Musical toy gives children an easy way to make music on their own simply by playing with blocks.</p> <p>Develops skills such as; shape and colour recognition, creativity, experimenting skills, music listening, manipulative skills, cause and effect.</p>	

<p>Musical Mat – Funky Footprints</p>	<p>Lights and sounds funky footprints. 10 footprints to walk up and 2 to land on lights, counting, funny sounds, and cheers.</p> <p>Great for getting interactive and inclusion fun.</p>	
<p>Sensory ball</p>	<p>Develops manipulative skills</p>	

<p>Sensory caterpillar</p>	<p>Encourages manipulative skills</p>	
<p>Switch Toy - Ariel Mermaid</p>	<p>A switch adapted toy that provides cause and effect when switch is activated by the child.</p> <p>Requires batteries and an additional switch for operation.</p>	

<p>Switch Toy - Bump and go car</p>	<p>A switch adapted toy that provides cause and effect when switch is activated by the child.</p> <p>Requires batteries and an additional switch for operation.</p>	
<p>Switch Toy - Buzz Lightyear</p>	<p>A switch adapted toy that provides cause and effect when switch is activated by the child.</p> <p>Requires batteries and an additional switch for operation.</p>	

<p>Switch Toy – Lion</p>	<p>A switch adapted toy that provides cause and effect when switch is activated by the child.</p> <p>Requires batteries and an additional switch for operation.</p>	 A plush lion toy with a large, shaggy orange-brown mane and a lighter orange-brown body. It has a white muzzle with a small black nose and two small black eyes. The lion is sitting upright on its haunches with its front paws resting on its knees. The background is plain white.
------------------------------	---	---



<p>Switch Toy - My Little Pony</p>	<p>A switch adapted toy that provides cause and effect when switch is activated by the child.</p> <p>Requires batteries and an additional switch for operation.</p>	
<p>Thomas the Tank Engine with Remote Control</p>	<p>Follow Me Thomas's controller has been switch adapted for him to chug along in a large 2m circle, with sound and music to capture the child's attention.</p>	



Switch Toy –
Nemo



A switch adapted toy that provides cause and effect when switch is activated by the child.


Requires batteries and an additional switch for operation.





<p>Switch Toy - Radio cassette player</p>	<p>A switch adapted toy that provides cause and effect when switch is activated by the child.</p> <p>Requires batteries and an additional switch for operation.</p>	 <p>The image shows a black portable radio cassette player. A custom switch panel is attached to the front, featuring a toggle switch labeled 'ADAPTED RADIO CASSETTE SWITCH' and 'PUSH TO ON / PULL TO OFF'. The panel also includes a 'VOLUME' knob and a 'PAUSE STOP' button.</p>
<p>Switch Toy - Rikki Rooster</p>	<p>A switch adapted toy that provides cause and effect when switch is activated by the child.</p> <p>Requires batteries and an additional switch for operation.</p>	 <p>The image shows a plush rooster toy with white and black feathers, a red comb, and orange feet. It is a stylized representation of a rooster.</p>

<p>Switch Toy – Robot</p>	<p>A switch adapted toy that provides cause and effect when switch is activated by the child.</p> <p>Requires batteries and an additional switch for operation.</p>	
<p>Switch Toy - Textured roller</p>	<p>This captivating device provides a very special somatosensory experience for children.</p> <p>Bumps on the switch offer tactile stimulation and help hand-eye coordination.</p> <p>Roller switch can be used as a switch-only when you spin the drum and the unit is activated - or as a switch-toy when you spin the drum to activate lights vibrations and music.</p> <p>Requires 2 AA batteries.</p>	


<p>Bongo Bright Beats switch adapted drum</p>	<p>Press your capability switch to hear up beat sounds, music & lights.</p> <p>Requires batteries and an additional switch for operation.</p>	
<p>Ring Around Bells</p>	<p>You can use this wonderful musical toy in four different ways: activate the attached yellow switch and the colourful bells twirl around, sounding the musical scale while the multi-coloured LEDs blink, or use a capability switch to activate the toy.</p> <p>Ring Around Bells is excellent for music therapy and music time. The bells are precision tuned and provide exceptionally clear sounds.</p> <p>Requires batteries for operation.</p>	


<p>Various weighted toys</p>	<p>Weighted stuffed animals are the ultimate sensory toy.</p> <p>The varied shapes, sizes and textures of our furry friends are not only appealing but they allow for weight to be distributed where it's needed most.</p> <p>Children of all ages and sizes to benefit from the calming effect of deep sensory, proprioceptive, and tactile sensory input.</p>	
------------------------------	---	---

9. Grab bars and suction rails


Grab bars and suction rails		
Stabilising horizontal Grab Bars and vertical Anchor Grab Bars	<p>The Stabilising Grab Bar encourages the child to bring the arms forward to stabilise, improving posture and increasing access to activities.</p> <p>It consists of a horizontal or vertical bar that attaches to the table top with suction cups, and can be used on the table or tray provided with a chair or standing frame.</p>	
Suction Rail	<p>Suction grab rails are lightweight and portable, featuring a durable commercial grade suction cup making them the strongest in the market.</p> <p>Easy to install and removes without any tools.</p>	



10. Positioning systems

Positioning systems		
Leckey Firefly Play Pak	<p>This lightweight bag contains an entire activity centre.</p> <p>Inside you will find rolls, wedges, and supports (designed by clinical experts) that stick together in dozens of different ways to let you and your child practice the early positions which form the building blocks for later abilities.</p> <p>Equipment measurements in bag: - H 20CM W 45CM L 45CM</p>	 A photograph of the Leckey Firefly Play Pak, a colorful, multi-colored bag (orange, yellow, green, red, blue) with a black strap and buckle, lying flat on a light surface. The bag is open, revealing various rolls and wedges inside.



<p>Wedge with detachable head rest</p>	<p>Allows patient to be positioned in prone - good positioning for eliminating the effect of gravity present when individual is upright.</p> <p>Allows patient to be positioned to stretch flexed hips and knees. Scoop shape decreases chance that child could roll off wedge.</p> <p>Headrest is useful for child with poor head and neck extension.</p> <p>Equipment measurements: - H 30 CM (at high end of wedge) W 43CM L 150CM</p>	 <p>The image shows a red and yellow wedge-shaped device, likely made of foam or plastic, used for patient positioning. It has a yellow top surface and a red side surface. Two blue straps are wrapped around the red surface, and a green strap is also visible. The device is placed on a wooden surface, and a metal shelf with various items is visible in the background.</p>
--	---	--


11. Ramps



Ramps		
Depac ramps	Portable folding ramps measuring 2 metres long	



<p>Easy access portable ramps</p>	<p>Easy access ramps are a lightweight yet portable ramp designed for doorways, sliding glass doors, and raised landings.</p> <p>One metre in length.</p>	
<p>Threshold ramp (small and large)</p>	<p>A small threshold ramp designed for use in doorways, sliding glass doors and raised landings.</p> <p>The Library has a variety of sizes stocked.</p>	



12. Accessories



Accessories		
Big keys keyboard	<p>Big keys keyboards for computer access provide users large, high-visibility keys that are easy to see.</p> <p>In addition, for users who lack fine motor movements, the larger keys provide a larger strike area that is easier to access compared to standard keyboards.</p>	
Campus SX and microphone boom with case	<p>Campus SX is a dedicated Phonak FM transmitter for educational applications when used in combination with Micro Boom.</p> <p>Professionals, teachers, educators and children benefit from maximum flexibility and user comfort.</p>	

<p>Phonak DigiMaster 5000 kit including Speaker with stand and Pass-around Mic</p>	<p>The Phonak's DigiMaster 5000 innovative loudspeaker array features no less than 12 individual high-quality loudspeakers, all housed in a robust aluminium frame.</p> <p>The DigiMaster 5000 array emits cylindrical sound waves at the most important speech frequencies, ensuring that the required sounds are carried effectively across a large area, while generating very little speech-hindering 'late reverberation' or echo. The DigiMaster 5000 is supplied with a sturdy floor stand</p> <p>The Pass-around Mic can be handed from one child to another, allowing the child who is wearing Roger Receivers to hear not only the teacher/educator but also other children's comments.</p>	 <p>The image displays the components of the Phonak DigiMaster 5000 kit. At the top is a tall, slender, black cylindrical loudspeaker array mounted on a black three-legged floor stand. To the left of the stand, a white mobile phone is connected to the system via a black cable. Below the phone, a coiled black cable with two small black circular components is visible. In the foreground, a white, handheld, cylindrical 'Pass-around Mic' with a green grille at the top is shown.</p>
--	---	---

<p>Platform steps</p>	<p>Platform to help assist children with posture when sitting and accessing benches.</p> <p>The Library stocks a variety of sizes.</p>	
<p>Steps with rails</p>	<p>Designed for assisted child access to change benches.</p> <p>Construction is lightweight and the handrails are adjustable in height.</p> <p>Equipment measurements: - H 73CM W 52CM L 70CM</p>	

<p>Vinyl wedges and shapes</p>	<p>The Library has a variety of sturdy and waterproof equipment that can promote inclusion and play activity.</p> <p>Types of equipment are wedges, various sized cylinders, different shaped animals including Loch Ness and Elephant, soft sided See Saw and double sided – 2 steps.</p>	
<p>Switch – Big Mac</p>	<p>Record any single message directly into the Big Mac and press its large activation surface for up to 2 minutes of playback.</p> <p>Students with visual disabilities and those who require a larger activation surface will have no problem connecting with BIG Mac communicator's large, 13cm (5-inch) activation surface as they connect with teachers, educators, children, and the world around them.</p>	

<p>Switch – Big Red</p>	<p>Large colourful surface makes this switch ideal for children with visual, cognitive, or physical requirements.</p>	
<p>Switch – Jelly Bean</p>	<p>Low profile circular switch on mounting plate for wheelchair tray or table top. Responsive across entire surface.</p> <p>Range of colours available. Switch can be mounted on small triangular plate, or two or four-place rectangular plate.</p> <p>Twist model enables switch tops to be interchanged with one of a range of colours.</p>	

<p>Switch - Smiley face</p>	<p>Switches enable children with severe disabilities to experience cause and effect, to play and have fun, and to interact with their environment.</p>	
<p>Various Auslan books, posters and flash cards</p>	<p>By teaching simple sign language to children from as young as seven months, we can help them to convey their emotions and their needs, as well as giving them a valuable head start in their social and intellectual development.</p>	

[illegible][illegible]